## 1. Allgemeine Informationen

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## 1. General Information

### Product description and Contents of Package

The Halcor® system serves to stabilise subcapital distal metatarsal-I-osteotomies when correcting hallux valgus and metatarsus primus varus deformities. The underlying technique corresponds with the Kramer technique or the minimally invasive technique according to Bösch.

The system consists of a correction plate with a T-profile and a screw thread as the head with 2 different offsets (4 or 6 mm). The T-shaped design of the fixation section achieves stability of the implant in metatarsal I. A fixed-angle spongy bone screw is used to fix the plate to the lateralised metatarsal head I.

 This product may only be used with the corresponding original system components and in accordance with the respective product-specific surgical technique and instructions for use.

### Information on reuse

Only unused products taken from the original packaging may be used. The HALCOR® plate and screw are not reusable.

### Properties of the material and biocompatibility

This product is made of a wrought titanium-aluminium 6-vanadium 4 alloy in accordance with ISO 5832-3. The entire implant surface has a matt gloss. Osteointegration is not intended.





### Warning / precautions

Warnings on the package must always be observed.

This product must only be used in its original condition. Modifications of any kind and mechanical processing are not permitted.

In general, the implants are sterilized with gamma radiation (at least 25 kGy). Implants that are delivered "sterile" must not be re-sterilized by the buyer.

### Key to the symbols used on the product label

	Do not reuse
	Do not re-sterilize
	Lot number
	Sterilization by irradiation
	Quantity
	Caution
	Order number
	Expiry date
	Manufacturer

## 2. Storage and handling

**Storage**  Sterile implants must always be stored unopened in their protective packaging. The implants are packed in such a way that they can be stored at normal room temperature and relative humidity (corresponding to normal "storage and working environment" in our latitudes; this means about 15°C and 35°C, depending on the time of year, and corresponding normal relative humidity) without risking any impairment of the packaging, the sterility, the product, etc. Do not expose sterile implants to sunlight without protection (i.e. without protection from light by the packaging carton).

**Handling**  Before using the sterile implant check the package for damage, which may have compromised product sterility.

If the original packaging is damaged the implant is no longer sterile and must not be implanted. When unpacking the implant, check that it corresponds with the description on the product label (part number and size). When removing the implant from the packaging, appropriate aseptic procedures must be observed by the surgical staff: The outer and, if applicable, the middle PE pouch can be removed by the non-aseptic surgical assistant; the innermost pouch must be removed by the aseptic surgical assistant. The implant must not come into contact with any objects that might damage its surface. Before insertion, every implant must be visually checked for damage.

## 3. Instructions for Use - Special notes on use

### Patient information

Every patient must be informed in detail about the surgery and about possible side effects prior to the operation. In particular, he should also be made aware that the safety and life-time of the implant can be influenced by his conduct (for instance by avoiding heavy physical activity and making sure that his body weight does not exceed the normal body weight for a person of his height). The surgeon should record all information given to the patient in writing.

### Preoperative planning

Preoperative planning delivers important information about the suitable component type and possible combinations. Additional implants should be kept at hand in case a different size is required or the selected implant cannot be used.

### Surgical technique

Implantation is essentially performed as described in the respective surgical technique. The surgeon must have sufficient training and experience in foot surgery, and must in particular be familiar with this system.

The following items must be observed in particular:

Do not use excessive and jerky force when using the HALCOR® plate as a lever to shift bone fragments.

The HALCOR® plate must always be tapped in with particularly sensitive, light taps on the hammer. Essentially, the excessive use of force when inserting or fixing the implant must be avoided due to the risk of chip fractures or bone avulsion.

### Indications

- Hallux valgus and metatarsus primus varus deformities with moderate to higher grade hallux deformities
- Deformities with a hallux rigidus component are particularly suitable
- Consequences of earlier surgery, e.g. relapsed hallux valgus

In the case of incongruent or fixed metatarso-phalangeal joints, soft tissue surgery is also performed.

### Relative indications

- Pseudarthrosis after subcapital metatarsal I osteotomies

### Contraindications

- Known intolerance to the implant material (Ti6Al4V) as well as individual components parts of it
- Acute or chronic infections, whether local or systemic.
- Severe muscular, neurological, or vascular disorders affecting the leg in question.
- Lack of bone substance or poor bone quality endangering stable fixation of the plate.

### Potential side-effects

The adverse effects listed below are among the most typical and frequent complications of implantation of the HALCOR® plate.

### Early complications

- Infection
- Haematoma
- Fracture
- Loosening of the implant, esp. the screw fixation
- Breakage of the implant
- Loss of correction
- Venous thrombosis and pulmonary embolism
- Intolerance reactions to the implant material

These side effects should be diagnosed and treated as soon as possible. They do not generally affect the outcome.

### Late complications

- Shift of position and loosening
- Breakage of the implant
- Pseudarthrosis
- Relapsed hallux valgus

The following conditions can have a negative effect on the outcome of surgery:

- Severe osteoporosis
- Severe deformations
- Local bone tumors
- Systemic diseases and metabolic disorders
- History of infections and falls
- Drug addiction or drug abuse
- Nicotine abuse
- Depressions
- Obesity
- Major physical activities associated with strong vibrations that expose the plate to impacts and/or excessive loads.

**Note** 

In the case of confirmed pseudarthrosis, rapid healing or revision with removal of the plate should be aimed for in order to prevent a fatigue fracture of the HALCOR® plate.

### Interactions with medicinal products

Interactions with medicinal products are unknown.

### Information about the Risks of Interfering Interaction

Undesired effects that harm the patient can occur when performing MRI/CT scans. Possible effects include artefacts, warning of the implant, induction of electrical currents, loosening of the implant. The manufacturer's instructions for use must be studied before the device is used. In case of doubt, reference implants must be tested for suitability in the respective MRI/CT instrument within the scope of an individual risk assessment. The patient must be informed about the risks.

### Postoperative treatment

- In accordance with the surgical technique
- Use recognized procedures and techniques in postoperative care and treatment.
- Document postoperative treatment in accordance with the hospital's standard procedures, e.g.:
  - surgical records
  - postoperative X-rays
  - regular checks / follow-up examinations.

## 4. Cleaning, disinfection and sterilization/re-sterilization of instruments

For information about the treatment of the instruments see the Cleaning and Treatment Instructions Lit. No. 860501. The manufacturer cannot be held liable for instrument cleaning and sterilisation processes carried out at the buyer's facilities.

### Important

Further information on the use of this product is available from your local sales office or on the respective website [www.implan-tec.at](http://www.implan-tec.at).

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