ImplanTec	<p>Prothesen für das ANA.NOVA® Hüftgelenk Prothèses pour articulation de hanche ANA.NOVA® Protesi d'anca ANA.NOVA® Hip joint prostheses for ANA.NOVA® Prótesis para la articulación de cadera ANA.NOVA®</p>
------------------	---

<i>Hersteller / Manufacturer</i>	<p>Phone: +43 2236 864 194 Fax: +43 2236 864 234 office@implan-tec.at www.implan-tec.at</p>
<i>ImplanTec GmbH</i>	
<i>Grenzgasse 38 a</i>	
<i>2340 Mödling</i>	

CE	0408	Lit. Nr. 841001
-----------	-------------	-----------------

1. Allgemeine Informationen

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1. General information

Product description and package contents

- The product, package contents, and the materials used are specified on the product labels. This product should only be used in accordance with its product-specific surgical technique.
- The batch and lot number(s) of the implants used must be recorded in the patient documentation. Labels are included in the packaging of the sterile implants for this purpose.




Information on reuse

- A prosthesis which has previously been implanted must not be reused.
- The implants are generally not intended for reuse. This rule, and others which may apply in exceptional circumstances, is specified by an entry on the product label.
- In general the implants are sterilized by means of gamma rays (at least 25 kGy) or ethylene oxide. Implants supplied "sterile" may not be re-sterilized by the buyer.

Warnings / precautions

- Warnings on the package must always be observed.
- This product must only be used in its original condition. It must not be altered or processed mechanically in any way, unless this is specifically required by the surgical technique.

Key to the symbols used on the product label

	Do not reuse
	Do not re-sterilize
	Lot number
	Sterilization by irradiation
	Sterilization by ethylene oxide
	Quantity
	Caution
	Order number
	Expiry date
	Cementless
	Manufacturer

2. Storage and Handling

Storage of the sterile implants

Implants must always be stored unopened in their protective packaging. The implants are packed in such a way that they can be stored at normal room temperature and relative humidity (corresponding to normal storage and working environment in our latitudes; this means about +15°C and +35°C, depending on the time of year, and corresponding normal relative humidity) without risking any impairment of the packaging, the sterility, the product, etc. Do not expose sterile implants to sunlight without protection (i.e. without protection from light by the packaging carton).

Handling of sterile implants

Before using the sterile implant check the package for damage, which may have compromised product sterility. If the original packaging is damaged the implant is no longer sterile and must not be implanted. When unpacking the implant, check that it corresponds with the description on the product label (part number and size).

When removing the implant from the packaging, appropriate aseptic procedures must be observed by the surgical staff: The outer and, if applicable, the middle PE pouch/blister can be removed by the non-aseptic surgical assistant; the innermost pouch/blister must be removed by the aseptic surgical assistant. The implant must not come into contact with any objects that might damage its surface. Before insertion, every implant must be visually checked for damage.

3. Instructions for use - special notes on use

Patient information

Prior to surgery, every patient must be informed in detail about the operation and possible adverse effects. In particular he/she should be advised that and how he/she may influence the safety and longevity of the implant (e.g. by avoiding excessive physical activity and by maintaining normal body weight). A general weight limitation is not intended.

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Preoperative planning

Preoperative planning provides essential information on the appropriate prostheses and likely component combinations. Additional implants should be kept at hand in case a different size is required or the selected implant cannot be used.

Surgical technique

The implantation must be carried out in accordance with the appropriate surgical technique, information on which can be obtained from ImplanTec GmbH. Familiarity with the recommended surgical technique for this system and its careful application are essential for optimal results.

Combination restriction

The implants and instruments manufactured and/or distributed by ImplanTec GmbH can be combined with each other.

The surgeon must always make sure that the individual implant components are compatible, and must take the general restrictions with regard to the combination of materials into account. Components made of stainless steel (ISO 5832-1, ISO 5832-9) must not be combined with components made of titanium (ISO 5832-3, ISO 5832-1).

Components made of forged stainless steel (ISO 5832-1) must not be combined with components made of cobalt chromium molybdenum alloy (ISO 5832-4) or forged components made of cobalt chromium molybdenum alloy (ISO 5832-12).

The components of the acetabular part of the joint (cup / inlay) and the femoral part of the joint (stem / head) should always be procured from a single supplier. All hip stems supplied by ImplanTec GmbH are made with a stem cone 12/14 manufactured to the specifications of CeramTec AG and are thus suitable for combination with the ceramic heads Biolox[®]-forte and Biolox[®]-delta (ISO 6474) S/M/L/XL.

Use of the femoral joint parts from one manufacturer within the meaning of the Medical Devices Act with the CeramTec AG heads may result in the possibility of combination with endoprostheses from other manufacturers. The combination with appropriate third-party products can be examined and compatibility assessed by ImplanTec GmbH on request.

Instructions for use

The hip stem cone and inner cone of the head must be clean and undamaged when assembled. The cone must be cleaned carefully before the head is fitted. The suitable head must then be fitted by hand and fixed on the cone using the fitting instrument and an appropriate stroke of the hammer.

The inner hip stem cone and outer cone of the inlay must be clean and undamaged when assembled. The inner cone must be cleaned carefully before the inlay is fitted.

The components of the hip endoprosthesis must be cleaned carefully before being assembled. The components must be clean and undamaged when assembled.

Indications

- Advanced hip joint destruction resulting from degenerative, posttraumatic, or rheumatoid arthritis
- Fracture or avascular necrosis of the femoral head
- Consequences of a previous surgery, e.g. osteosynthesis, joint reconstruction, arthrodesis, hemi-arthroplasty, or total hip replacement

Contraindications

- Acute or chronic infections, whether local or systemic
- Severe muscular, neurological, or vascular disorders affecting the leg in question
- Lack of bone substance or poor bone quality endangering stable fixation of the prosthesis
- Any concomitant illness which could compromise functioning of the prosthesis
- Revision with extensive bone defects

Possible adverse effects

The adverse effects listed below are among the most typical and frequent complications of hip arthroplasties.

Early complications:

- Local intolerance and loosening of the prosthesis
- Toxic systemic side effects and
- Allergies to material components
- Venous thrombosis and pulmonary embolism
- Cardiovascular and pulmonary (e.g. fat embolism) dysfunctions
- Haematoma

These side effects should be diagnosed and treated as soon as possible. They do not generally affect the outcome.

Late complications:

- Migration and loosening of the prosthesis
- Abrasion of the prosthesis surface and development of osteolysis due to a foreign body reaction.

The following conditions can impair the outcome of the operation:

- Severe osteoporosis
- Severe deformities, congenital hip dislocation
- Local bone tumours
- Systemic and metabolic disorders
- History of infections and falls
- Drug addiction or abuse

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- Nicotine abuse

- Depressions
- Obesity

- Heavy physical activity or activities associated with severe jarring in which the prosthesis could be subject to impacts and/or excessive stresses (e.g. heavy physical labour, marathon running, etc.).

Avoidance/minimization of side effects

If the implantation of this hip prosthesis is regarded as the best solution for the patient while some of the above conditions apply, it is very important to inform the patient of the expected effects of these circumstances on the success of the operation. It is also recommended to advise the patient on how to reduce the impact of adverse circumstances. The surgeon should record all information given to the patient in writing.

Postoperative treatment

Use recognized procedures and techniques in postoperative care and treatment. Document postoperative treatment in accordance with the hospital's standard procedures, e.g.:

- Surgical records
- Postoperative X-rays
- Patient's card
- Regular checks / follow-up examinations

Cleaning, disinfection and sterilization/re-sterilization of instruments

The following recommendations are for information only. The manufacturer accepts no liability for the cleaning and sterilizing processes of instruments performed in the purchaser's facilities.

Instructions for hygiene treatment of instruments

The instruments are supplied "non-sterile" and must be cleaned, disinfected and sterilized using validated processes before use.

The following instructions are based on validation trials conducted in the course of preparing the medical device for reuse. The person performing the hygiene treatment is responsible for the treatment actually provided achieving the desired results using the equipment and people available in the hygiene treatment facility.

For precise details about preparation, preliminary manual cleaning, checking & maintenance, packing and storage, please see our hygiene treatment instructions, ref. no. 860501-D!

Cleaning and disinfecting by machine

Cleaning and disinfecting by machine should be performed by a process validated by the relevant regulations in a washer-disinfector complying with ÖNORM EN ISO 15883-1 and -2, preferably using a (mild) alkaline cleaning agent. The following program was used for prevalidation:

Phase	H ₂ O type	Temp. [°C]	Time [min]	Dosing temp. [°C]	Dosing [ml/l]
Initial rinse	purified	cold	≥ 2	-	-
Wash	purified	55+/- 5	≥ 7	35	4 - 10
Rinse I	purified	cold	≥ 1	-	-
Rinse II	purified	cold	≥ 1	-	-
Disinfect	purified	93 +/- 3	≥ 5	-	-
Dry	-	110	15	-	-

Sterilizing
Steam sterilizing with a process validated to ÖNORM EN ISO 17665/ONR 112069 in a steam sterilizer with a fractionated vacuum process as per ÖNORM EN 285. Prevalidation was conducted with the following program specifications:

Cycle type	Temperature
------------	-------------