1. Allgemeine Informationen

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1. General Information

Product description and package contents

The CCG system consists of CCG and CCG-GF bands as well as CCG stabilisers. CCG bands and CCG stabilisers are patented products. The CCG band is formed into a loop and compressed using a special cerclage. The CCG stabiliser is a curved plate with dentiform protrusions on both sides which penetrate into the bone, thus providing a movement-free bone support. CCG stabilisers must only be used together with CCG bands tied around them. CCG bands and CCG stabilisers are made of pure titanium according to ISO 9001/46001 and ASTM BF 67 Gr. 4. The titanium surface is roughened, enabling osseointegration. The batch or lot number(s) of the implants used must be documented in the patient records. Labels are included in the packaging of the implants for this purpose.

Information on reuse

A previously used band which has had to be removed must not be reused because of the deformation which has occurred. In this case, the cerclage must be repeated using a new band.

Warnings and precautions

These products may only be used in their original condition or with modifications expressly described in the surgical technique. For the successful implantation of CCG compression bands and CCG stabilisers it is of course a requirement that the surgeon be fully conversant with the methods of osteosynthesis when treating bone fractures and the problems associated with revision operations of both cemented and cementless hip arthroplasties.

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	Do not reuse
	Do not re-sterilize
	Lot number
	Sterilization by irradiation
	Quantity
	Caution
	Order number
	Expiry date
	Manufacturer

2. Storage and handling

Storage

Sterile implants must always be stored unopened in their protective packaging. The implant packaging is designed so that it can be stored at normal room temperature and relative humidity (corresponding to normal storage and working environment in our latitudes; this means about +15°C and +35°C, depending on the time of year, and corresponding normal relative humidity) without risking any impairment of the packaging, the sterility, the product, etc. Do not expose sterile implants to sunlight without protection (i.e. without protection from light by the packaging carton).

Handling

CCG Bands

Before use, the packaging of the sterile implant must be examined for damage, as this can affect the product sterility. When removing the implant from the packaging, appropriate aseptic procedures must be observed by the surgical staff: The outer PE pouch is opened by the non-aseptic surgical assistant; the innermost pouch is opened by the aseptic surgical assistant.

CCG stabiliser

The outer PE pouch is opened by the non-aseptic surgical assistant; the innermost pouch is opened by the aseptic surgical assistant.

CCG bands and CCG stabilisers

The implant must not come into contact with any objects that might damage its surface. Before insertion, every implant must be visually checked for damage.

3. Instructions for Use - Special notes on use

Patient information

Prior to surgery, every patient must be informed in detail about the operation and possible adverse effects.

A. Surgical technique CCG bands

1. Reposition the bone fragments to be fixed using reposition forceps. Repositioning must not be carried out with the CCG band.
2. Using a curved bone rasp, expose the bone subperiostally at the cerclage position.
3. Insert the CCG loopawl from above so that it maintains close contact with the bone. If using the CCG-GF: Slide the moveable spike into the appropriate position.
4. Hook in the CCG band so that the clasp is on the outside and pull through.
5. Pull CCG compression band through the clasp and tighten.
6. Bring the carriage of the CCG cerclage tightener all the way to the front by turning to the left, then tighten the band and fixate it with the lever.
7. The compression is carefully carried out by turning to the right with thumb and index finger only. The intended compression has been achieved when the band cannot be moved.

Caution: Using excessive force can result in undesirable over-compression! When the gap of the clasp begins to open, the permissible compression has already been slightly exceeded. The compression must not be continued at this point.

8. Bend the cerclage tightener with the band slowly up to 90°.

WARNING

During bending with the cerclage tightener, a combination of tensile and bending forces affect the CCG band. Therefore the bending may only be performed up to 90° so that the band will not slip back. Quick bending over 90° will increase the tensile and bending forces, which increases the strain, and the CCG band may break. Once the CCG cerclage tightener has been removed, the bending is completed by the cerclage fixator according to the instructions, using only bending, but not tensile forces. This eliminates the risk of breakage. This handling instruction is very important and must be followed.

9. Loosen the lever and remove the cerclage tightener.
10. Insert the CCG fixator into the CCG band, bend the band, and shorten it. Use of the CCG fixators eliminates the risk of band breakage during bending, and it prevents an incorrect shortening of the band. The wire cutter is not a substitute for the CCG fixator.
11. Set the fixator with one end onto the end of the band and impact it firmly. For application with wire cutter and pestle: Instead of the CCG fixator, the CCG band is cut off by a wire cutter 1 cm above the clasp. Then the protruding end of the band is bent by a pestle and impacted onto the clasp.

B. Surgical technique CCG stabiliser

1. Only use the CCG stabiliser in combination with CCG cerclage bands.
2. Choose the correct length of stabiliser according to the cortical bone to be supported or the shape of the fracture line.
3. The stabiliser is applied longitudinally, axially to the bone, and secured using reduction forceps. Bone reduction should be carried out using the reduction forceps, and not by the compressive effect of the CCG bands.
4. The number of stabilisers used will depend on the cortical defect to be bridged. The positioning and spacing of the stabilisers will have to be adapted to the intra-operative situation.
5. When using several CCG stabilisers, these should be placed parallel to one another. Contact between the stabilisers must be avoided in all cases.
6. Furthermore, the stabiliser should be positioned in such a way that direct contact with other metal implants is excluded.
7. Use a tappet to tap the stabilisers protrusions evenly into the bone.
8. Secure the stabiliser into place with CCG bands.
9. When using two CCG cerclage bands, position one band proximally, the other distally to the reposition forceps. When using more than two CCG cerclage bands or longer stabilisers, the use of an additional reposition forceps is recommended.
10. The CCG cerclage bands must be placed perpendicularly to the stabiliser between two grooves to ensure optimum connection between CCG cerclage band, stabiliser and bone.

The CCG system is suitable for particular forms of osteosynthesis, particularly for cerclage of the greater trochanter and femur associated with revision to hip arthroplasties. In contrast to cerclage wires, the CCG compression band allows a functional compression to be achieved. It is assumed that an implant revision is already indicated in this connection.

The CCG system fulfils a supporting function. The following indications apply:

- Reattachment of the greater trochanter
- Bone support in case of a damaged proximal femur (e.g. by bone cement, loosened stem, or polyethylene wear particles), with a possible cancellous bone graft in the proximal femur area.
- Cerclage of bone fissures.
- Fracture protection during rasping of a new prosthetic canal and prior to implantation of a new stem.
- Closing of cortical fenestration.
- Fractures in the area of the femoral prosthesis stem.
- To improve bone-implant contact.

In general, CCG and CCG-GF have the same function. In the case of conical shape of the bone, primarily above the lesser trochanter, the CCG-GF is indicated rather than the CCG. CCG stabilisers are specially indicated in fractures and shaft fissures in which the use of CCG bands alone would not grant sufficient stability.

Contraindications

An implant of suitable design and length is a prerequisite for a successful use of the CCG system in hip revisions. Complications such as unstable implantations, implant subsidence, stem and prosthesis breakage cannot be prevented entirely by the CCG system. Therefore CCG is not indicated in such cases. The CCG system is not indicated in cases with large defects of the cortical bone and hence unavoidable contact with the implant which would lead to metal wear.

Complications

CCG band may break as a result of improper handling. Please read the instructions. Complications which result from incorrect indications, surgical technique or aseptic conditions remain the responsibility of the operating surgeon and cannot be attributed to either the manufacturer or the supplier of CCG.

Avoidance and minimization of side effects

A tolerance test should be performed before the operation in the case of any uncertainties about the patient's tolerance of the implant material. The appropriate allergy test patches and test instructions are available from your local sales office.

Postoperative treatment

In accordance with the surgical technique. Use recognized procedures and techniques in postoperative care and treatment.

Document postoperative treatment in accordance with the hospital's standard procedures, e.g.:

- surgical records
- postoperative X-rays
- patient's card
- regular checks and follow-up examinations.

4. Cleaning, disinfection and sterilization/re-sterilization of instruments

The following recommendations are for information only. The manufacturer accepts no liability for the cleaning and sterilizing processes of implants and instruments performed in the purchaser's facilities. In general the implants are sterilized by means of gamma rays (at least 25 kGy). Implants supplied "sterile" may not be re-sterilized by the buyer.

Instructions for hygiene treatment

The CCG instruments are supplied "non-sterile" and must be cleaned, disinfected and sterilized using validated processes before use. The following instructions are based on validation trials conducted in the course of preparing the medical device for reuse. The person performing the hygiene treatment is responsible for the treatment actually provided achieving the desired results using the equipment and people available in the hygiene treatment facility. For precise details about preparation, preliminary manual cleaning, checking & maintenance, packing and storage, please see our hygiene treatment instructions, ref. no. 860520 D.

Cleaning and disinfecting by machine

Cleaning and disinfecting by machine should be performed by a process validated by the relevant regulations in a washer-disinfector complying with ÖNORM EN ISO 15883-1 and -2, preferably using a (mild) alkaline cleaning agent. The following program was used for prevalidation:

Phase	H ₂ O-Art	+ H ₂ O (L)	°C	min
Initial rinse	Cold	0	-	2
Wash	Demineralized	+ 3,0	55	7
Neutralize	Hot	0	W	2
Rinse	Hot	0	W	1
Disinfect	Demineralized	0	93	5
Dry	-	-	110	15

Sterilizing

Steam sterilizing with a process validated to ÖNORM EN ISO 17665/ONR 112069 in a steam sterilizer with a fractionated vacuum process as per ÖNORM EN 285. Prevalidation was conducted with the following program specifications:

Cycle type	Temperature	Pressure	Sterilization time	Drying time
Fract. vacuum	134 °C	3050 mbar abs.	5 minutes	25 minutes

Other sterilizing methods/cycles can be used, but these should be validated in accordance with the relevant regulations before being used routinely.

Important

Further information on the use of this product is available from your local sales office. The medical author of the CCG is available to provide instruction to interested surgeons.

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